

Clinisept+

## Cautions:

Always follow the instructions carefully. For external use only. Avoid contact with eyes. If contact with eyes occurs, rinse with running water. Do not ingest. If ingested, seek medical advice as a precautionary measure. Clinisept+ has been dermatologically tested and proven to be non-irritant. Cease usage if irritation occurs. Avoid direct contact with clothing as may discolour fabric.

This product requires no specific disposal procedure. Use within 3 months of opening. Clinisept+ may be decanted into a gallipot for immediate use but do not decant into another container.

Avoid extremes of temperature and keep out of direct sunlight.

All components fully recyclable - please recycle responsibly.

## Ingredients:

Water (Aqua), Phosphoric Acid, Sodium Hypochlorite.

Contains no chlorhexidene, alcohol or fragrance and is colourant free.

Visit [www.cliniseptplus.com](http://www.cliniseptplus.com) for more information.



CLINICAL HEALTH  
TECHNOLOGIES

Made in the UK

Clinical Health Technologies Ltd, Hinckley, LE10 3BP

+44(0)1455 618866

[www.cliniseptplus.com](http://www.cliniseptplus.com)

Clinisept+

## Prep & Procedure

Fast Acting Skin Cleanser  
& Antimicrobial Solution

Dermatologically Tested

User Guide

## Clinisept+ Prep & Procedure

Please read this leaflet carefully to ensure that you use **Clinisept+ Prep & Procedure** correctly to enable it to deliver the best possible results.

**Clinisept+ Prep & Procedure** is a gentle, but highly effective skin cleanser that **cleans, calms and cares** for the skin. It contains the latest hypochlorous technology to ensure rapid antimicrobial protection, but without stinging or irritation. Because of its unique oxidising method of action it is not susceptible to antimicrobial resistance.

**Clinisept+ Prep & Procedure** is unlike any other cleanser and represents a significant development in skin care. It provides optimal conditions for skin recovery:

- Cleansing antimicrobial
- Skin neutral pH
- Non-irritant
- Non-cytotoxic
- Dermatologically tested
- No Chlorhexidine or Alcohol

**Clinisept+ Prep & Procedure** is suitable for use on normal and sensitive skin.

For professional use only.

## Using Clinisept+ Prep & Procedure:

Apply **Clinisept+ Prep & Procedure** by wetting a sterile gauze or cleansing pad with Clinisept+ and then gently wiping over the area using a circular motion. The volume of Clinisept+ applied should be adjusted to suit the size of the area to be cleansed. In all instances, ensure the whole area is thoroughly wetted. Sterile gauzes and pads should be replaced and discarded as necessary.

Once applied, **Clinisept+ Prep & Procedure** leaves no residue on the skin and can be left to dry naturally. Do not dry the skin after application unless with a gauze, dressing pad or clean towel.

Clinisept+ is rapid in action and requires a wet contact time of only 15 seconds.



## Cleansing before, during and after a procedure:

A skin neutral pH and non-cytotoxic chemistry allows Clinisept+ to be used before, during and after a treatment, thereby enabling high levels of hygiene to be maintained throughout a procedure. Use fresh sterile gauzes or cleansing pads as required and discard used pads.

**Clinisept+ Prep & Procedure** provides a cooling and calming sensation and is very effective in reducing redness. For maximum care of recovering or irritated skin, (ie. immediately following a treatment or procedure), thoroughly soak a gauze mask/square or cleansing pad with Clinisept+ and leave it placed on the affected area for five to ten minutes. Repeat as necessary.

